#### English

The ResMed Oxygen Connector Port (26960) is an accessory validated for use with up to 15 Litres per minute of supplemental oxygen with the following products: VPAP™ III ST-A, VPAP™ III ST-A with QuickNav, VPAP™ Auto 25, VPAP™ ST, VPAP™ S, VPAP™ IV and VPAP™ IV ST.



The Oxygen Connector Port is designed to allow the addition of supplemental oxygen at the flow generator output while reducing the amount of oxygen backflow into the flow generator compared with standard “T-piece” oxygen fittings.

#### Fitting the Oxygen Port

1. Firmly connect the Oxygen Connector Port (A) directly to the air outlet (B) of the flow generator or humidifier as shown.
2. Connect the air tubing to the end of the Oxygen Connector Port (C) as shown. Ensure that the tube is connected up to the indicated line (D).
3. Connect the oxygen supply tubing to the oxygen inlet port (E) as shown.

#### Notes:

- The arrows on the Oxygen Connector Port should be pointing along the direction of the airflow—away from the flow generator.
- At a fixed rate of supplemental oxygen flow, the inhaled oxygen concentration varies, depending on the pressure settings, patient breathing pattern, mask, point of application and leak rate.
- When oxygen is added, always ensure that you verify the correct operation of triggering and cycling, and activation of mask alarm when enabled.

#### Cleaning

- Wash the Oxygen Connector Port in warm water using mild detergent.
- Rinse thoroughly and allow to dry.

#### Disposal

The Oxygen Connector Port does not contain any hazardous substances and may be disposed of with normal household refuse.

#### Symbols which appear on the Product

Arrow indicates direction of airflow.

#### WARNINGS

- Always use the ResMed Oxygen Connector Port when adding supplemental oxygen at the flow generator.
- Oxygen flow must be turned off when the flow generator is not operating, so that unused oxygen does not accumulate within the flow generator enclosure and create a risk of fire.
- Oxygen supports combustion. Oxygen must not be used while smoking or in the presence of an open flame.
- Do not use bleach, chlorine, alcohol, or aromatic-based solutions, moisturising or antibacterial soaps or scented oils to clean the Oxygen Connector Port. These solutions may react with the plastic and reduce the life of the product.

#### CAUTION (USA ONLY)

Federal law restricts these devices to sale by or on the order of a physician.

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© 2012 ResMed Ltd. 268322/1 2012-01
S8 and VPAP are trademarks of ResMed Ltd. AutoSet, SmartStart and VPAP are registered in U.S. Patent and Trademark Office