

User Guide

A

B

English

Please ensure that you read any specific references in this user guide that apply to your product type. Use the following list to identify your product type:

- APAP—S8™ AutoSet™ devices
- CPAP—all other S8 devices (except S8 AutoSet)
- Bilevel—VPAP™ IV and VPAP™ IV ST devices.

Indications for Use

Your device is intended for home and hospital use.

CPAP

Your continuous positive airway pressure (CPAP) device is indicated for the treatment of obstructive sleep apnoea (OSA) in adult patients.

APAP

Your self-adjusting sleep apnoea device is indicated for the treatment of obstructive sleep apnoea (OSA) in adult patients. The self-adjusting sleep apnoea device has two treatment modes: AutoSet and fixed-pressure CPAP.

Bilevel

Your bilevel ventilator device is intended to provide noninvasive ventilation for patients with respiratory insufficiency or obstructive sleep apnoea (OSA).

Contraindications

Positive airway pressure may be contraindicated in some patients with the following conditions:

- pneumothorax or pneumomediastinum
- pathologically low blood pressure, particularly if associated with intravascular volume depletion
- cerebrospinal fluid leak, recent cranial surgery or trauma
- severe bullous lung disease
- dehydration.

Adverse Effects

You should report unusual chest pain, severe headache or increased breathlessness to your prescribing physician. An acute upper respiratory tract infection may require temporary discontinuation of treatment.

The following side effects may arise during the course of therapy with the flow generator:

- drying of the nose, mouth or throat
- bloating
- ear or sinus discomfort
- eye irritation
- skin rashes
- nosebleed.

Masks and Humidifiers

Recommended masks and humidifiers are available from your place of purchase. For information on using your mask or humidifier, refer to the manual supplied with your mask or humidifier.

Setting up your Flow Generator

Refer to illustration A on page 1.

1. Connect the power cord to the socket at the rear of the flow generator.
2. Plug the other end into a power outlet.
3. Connect the air tubing to the air outlet of the flow generator.
4. Connect the assembled mask system to the free end of the air tubing.

WARNING

- Make sure the power cord and plug are in good condition and the equipment is not damaged.
- Only ResMed air tubing should be used with the device. A different type of air tubing may alter the pressure you actually receive and reduce the effectiveness of your treatment.

- Blocking the tubing and/or air inlet of the device while in operation could lead to overheating of the device.
- Always use the ResMed Oxygen Connector Port for a VPAP IV or VPAP IV ST device when adding supplemental oxygen at the flow generator outlet.

CAUTION

- Be careful to place the device where it cannot be bumped, and where no one will trip over the power cord.
- If you put the device on the floor, make sure the area is free from dust and clear of bedding, clothes or other objects that could block the air inlet.
- Make sure the area around the device is dry and clean.

Flow Generator Control Panel

Refer to illustration B on page 1.

The control panel has the following items: (1) Start/Stop; (2) LCD screen; (3) Up and Down keys; (4) Left key; (5) Right key.

Starting Treatment

1. Make sure the power is on. If fitted, turn the main power switch at the back of the unit to on (I).
2. Lie down and arrange the air tubing so that it is free to move if you turn in your sleep.
3. To start treatment, press , unless your device has been set up to start delivering therapy automatically (SmartStart™/Stop). If so, your device will start when you breathe into your mask.
4. Fit your mask as described in the mask user instructions.


WARNING

A mask should not be used unless the flow generator is turned on and operating properly.

CAUTION

Do not leave long lengths of air tubing around the top of your bed. They could twist around your head or neck while you are sleeping.

Stopping Treatment

To stop treatment at any time, remove your mask and press , or, if your device has SmartStart/Stop enabled, simply remove your mask and treatment will end.

Using the Menus

Key Functions

Up and Down keys: Scrolls through menus and setting options.

Left key: Enters menus; changes options; applies settings.

Right key: Exits menu; cancels operation.

Menu Settings

Some menu items may not be available on your model.

Ramp: Period during which the flow generator increases from a low pressure to the prescribed treatment pressure.

Setting: Period during which the flow generator maintains a minimum pressure before adjusting pressure in response to respiratory events.

EPR Level: Expiratory pressure relief can be used to ease any discomfort when breathing out during treatment. Options: OFF (no EPR) – 3 (highest EPR).

Altitude: Selects the altitude range.

Mask: Selects your mask type.

Tube Length: Selects the length of air tubing.

Humidifier: Selects the type of humidifier used with your flow generator.

SmartStart/Stop: Turns the SmartStart/Stop function on or off.

Leak Alert/Mask Alarm: If enabled, sounds an alarm when a high mask leak is detected.

Language: Selects the LCD display language.

Cleaning and Maintenance

Daily:

Disconnect the air tubing and hang it in a clean, dry place until next use.

Weekly:

1. Remove the air tubing and wash it in warm water and mild detergent. Rinse thoroughly, hang and allow to dry.
2. Before next use, reassemble the mask and headgear.
3. Reconnect the air tubing.

Monthly:

1. Clean the exterior of the flow generator with a damp cloth.
 2. Check the air filter for holes and blockage by dirt. Replace it every six months, or more often in a dusty environment.
- For information on cleaning your mask or humidifier, refer to the manual supplied with your mask or humidifier.

WARNING

Beware of electric shock. Do not immerse the flow generator or power cord in water. Always unplug the flow generator before cleaning and make sure it is dry before reconnecting.

CAUTION

- Do not hang the air tubing in direct sunlight as it may harden over time and eventually crack.
- Do not use bleach, chlorine, alcohol, or aromatic-based solutions, moisturising or antibacterial soaps or scented oils to clean the air tubing or the flow generator. These solutions may cause hardening and reduce the life of the product.

Changing the Air Filter

1. Disconnect the power cord and air tubing.
2. Remove the air filter from the air inlet at the back of the flow generator.
3. Put the new air filter into the air inlet with the blue side facing out from the flow generator.

CAUTION

Do not wash the air filter.

Technical Specifications

Operating pressure range	4 to 20 cm H ₂ O (CPAP/APAP) 2 to 25 cm H ₂ O (Bilevel)
Power supply	Input range 100–240V, 50/60Hz S8: 2.5A < 140VA (110W) (maximum power consumption). VPAP: 40VA (typical power consumption), < 100VA (maximum power consumption).
Operating temperature	+5°C to +35°C
Operating humidity	10–95% non-condensing
Storage and Transport Temperature	-20°C to +60°C
Storage and Transport Humidity	10–95% non-condensing
Atmospheric pressure range	1,060 hPa (sea level) to 680 hPa (2,591 m)

Additional Information

Additional technical specifications and information can be found in the device-specific user guide on www.resmed.com on the **Products** page under **Service & Support**.

Electromagnetic Compatibility

Product complies with all applicable electromagnetic compatibility requirements (EMC) according to IEC60601-1-2, for residential, commercial and light industry environments. The electromagnetic compatibility tables for these ResMed devices can be found on www.resmed.com on the **Products** page under **Service & Support**. Click on the PDF file for your device.

Aircraft Use

Please consult ResMed and the medical services department of your carrier if you intend to use your flow generator on an aircraft.

The manufacturer reserves the right to change these specifications without notice.

Symbols that may Appear on the Product

-  Attention, consult accompanying documents
-  Type CF equipment
-  Class II equipment
-  Drip proof
-  Name and address of manufacturer
-  Dangerous voltage

 **Environmental Information** This device should be disposed of separately, not as unsorted municipal waste. To dispose of your device, you should use appropriate collection, reuse and recycling systems available in your region. The use of these collection, reuse and recycling systems is designed to reduce pressure on natural resources and prevent hazardous substances from damaging the environment. If you need information on these disposal systems, please contact your local waste administration. The crossed-bin symbol invites you to use these disposal systems. If you require information on collection and disposal of your ResMed device please contact your ResMed office, local distributor or go to www.resmed.com/environment.

General Warnings and Cautions

Warnings

- Read the entire guide before using the flow generator.
- The flow generator is NOT a life support ventilator and may stop operating with power failure or in the unlikely event of certain fault conditions.
- Advice contained in this guide should not supersede instructions given by the prescribing physician.
- A patient should not connect a device to the data communication port unless instructed to do so by their health care provider or physician. Only ResMed products are designed to be connected to the data communication port. Connecting other devices could result in injury, or damage to the flow generator. (Data communication ports are not available on all devices.)
- The flow generator should be used with masks (and connectors)* recommended by ResMed, or by a physician or respiratory therapist. A mask should not be used unless the flow generator is turned on and operating properly. The vent hole or holes associated with the mask should never be blocked.

Explanation: The flow generator is intended to be used with special masks (or connectors) which have vent holes to allow continuous flow of air out of the mask. When the device is turned on and functioning properly, new air from the device flushes the exhaled air out through the mask vent holes. However, when the device is not operating, insufficient fresh air will be provided through the mask, and the exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can, in some circumstances, lead to suffocation. This applies to most models of positive airway pressure devices.

* Ports may be incorporated into the mask or in connectors that are near the mask.

- In the event of power failure or machine malfunction, remove the mask.
- In certain fault conditions, pressures up to 30 cm H₂O (CPAP/APAP) or 40 cm H₂O (Bilevel) are possible.
- If oxygen is used with this device, the oxygen flow must be turned off when the device is not operating. If oxygen has been left on, turn off the device, then wait 30 minutes before turning on the device again.

Explanation: When the flow generator is not in operation and the oxygen flow is left on, oxygen delivered into the air delivery

tubing may accumulate within the flow generator enclosure and create a risk of fire. This applies to most types of positive airway pressure devices.

- Oxygen supports combustion. Oxygen should not be used while you are smoking or in the presence of an open flame.
- Always ensure airflow is being generated by the device before the oxygen supply is turned on.
- Always turn the oxygen supply off before stopping the airflow from the device.

Note: At a fixed rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on where the oxygen is introduced, the pressure settings, patient breathing pattern, mask selection, and leak rate.

- Do not use the flow generator if there are obvious external defects, unexplained changes in performance or unusual noises.
- Do not open the flow generator case. There are no user serviceable parts inside. Repairs and internal servicing should only be performed by an authorised service agent.
- Explosion hazard—do not use in the vicinity of flammable anaesthetics.
- The flow generator should not be used with anaesthetised patients.
- The device should not be connected to both AC and DC power sources simultaneously.
- In the clinical environment any PC that is used with the flow generator must be at least 1.5 m away from, or at least 2.5 m above the patient. It must also comply with IEC 60950 or equivalent.

Cautions

- At low pressures, the flow through the exhalation ports of your mask may not clear all exhaled gas from the tubing. Some rebreathing may occur.
- The air temperature for breathing produced by this device can be as much as 6°C higher than the temperature of the room. Caution should be exercised if the room temperature is warmer than 32°C.
- (S8 series only) When AC mains power (100–240V AC) is not available, always use a ResMed DC-12 converter. (The DC-12 converter is available as an optional accessory. It is not supplied with all models.)
- Do not remove any attached accessories while power is connected to your device.

Note: The above are general warnings and cautions. Specific warnings, cautions, and notes appear with the relevant instructions in this user guide.

Servicing

The flow generator should be inspected by an authorised ResMed service centre five years from the date of manufacture. Before this, the device is intended to provide safe and reliable operation if it is operated and maintained according to the instructions provided by ResMed. Warranty details are provided with the device at the time of original supply. As with all electrical devices, if any irregularity becomes apparent, you should have the device inspected by an authorised ResMed service centre.

Limited Warranty

ResMed warrants that your ResMed flow generator shall be free from defects in material and workmanship for a period of two years from the date of purchase by the initial consumer. This warranty is not transferable.

If the product fails under conditions of normal use, ResMed will repair or replace, at its option, the defective product or any of its components. This Limited Warranty does not cover:

- a) any damage caused as a result of improper use, abuse, modification or alteration of the product;
- b) repairs carried out by any service organisation that has not been expressly authorised by ResMed to perform such repairs;

- c) any damage or contamination due to cigarette, pipe, cigar or other smoke;
- d) any damage caused by water being spilled on or into the product.

Warranty is void on product sold, or resold, outside the region of original purchase. Warranty claims on defective product must be made by the initial consumer at the point of purchase.

This warranty is in lieu of all other express or implied warranties, including any implied warranty of merchantability or fitness for a particular purpose. Some regions or states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

ResMed shall not be responsible for any incidental or consequential damages claimed to have occurred as a result of the sale, installation or use of any ResMed product. Some regions or states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation may not apply to you. This warranty gives you specific legal rights, and you may also have other rights which vary from region to region. For further information on your warranty rights, contact your local ResMed dealer or ResMed office.

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ResMed Ltd (Manufacturer) 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia • **ResMed Corp** (US Designated Agent) 14040 Danielson Street Poway CA 92064-6857 USA • **ResMed (UK) Ltd** (EU Authorised Representative) 96 Milton Park Abingdon Oxfordshire OX14 4RY UK • **ResMed Offices** Australia, Austria, Belgium, Brazil, China, Finland, France, Germany, Greece, Hong Kong, India, Ireland, Italy, Japan, Malaysia, Mexico, Netherlands, New Zealand, Norway, Portugal, Singapore, South Africa, Spain, Sweden, Switzerland, UK, USA (see www.resmed.com for contact details).

Protected by patents and design registrations.

AutoSet, S8, SmartStart and VPAP are trademarks of ResMed Ltd. AutoSet, SmartStart and VPAP are registered in U.S. Patent and Trademark Office

© 2008 ResMed Ltd. 228262/2 08 07