
Chin Restraint Nasal Mask Accessory

USER'S GUIDE

INTENDED USE

The intent of the Chin Restraint is to restrict the mouth from opening during Positive Airway Pressure (PAP) treatment and minimise air loss through the mouth. It can be used in conjunction with a wide range of mask and pillow systems.

If possible, have the Chin Restraint initially fitted by staff at your sleep disorders center or by your equipment supplier.

FITTING THE CHIN RESTRAINT

The Chin Restraint can be fitted either before or after fitting the mask.

- 1 Fit the Chin Restraint over your chin so that the edge with the stitched 'V' mark is closest to your bottom lip. Bring the straps up past your cheeks and fasten the Velcro® at the top of your head.
- 2 If necessary, adjust the straps to provide sufficient support to your chin.

If you require a more secure fit, thread the Chin Restraint under your headgear strap as follows:

- 1 Fit the mask and headgear, then fit the Chin Restraint over your chin.
- 2 Pass the Chin Restraint strap on your right-hand side (without the Velcro tab) under the headgear strap at the top center of your head.
- 3 Then pass the other strap (with the Velcro Tab) across the top of both the headgear and Chin Restraint straps. Fasten the Velcro Tab to lock both the Chin Restraint and headgear in position.

Note: *If you experience any difficulty with the fitting of the Chin Restraint, consult your physician or equipment supplier.*

CAUTION

- Do not position the Chin Restraint so that it obstructs your mouth.
- Do not position the Chin Restraint completely beneath the chin or across your upper throat.
- Do not overtighten the Chin Restraint.
- If symptoms of a dry/sore mouth/throat occur, consult your physician or sleep therapist.
- If any visible deterioration (eg, wear, fraying, excessive stretching, tears) of the Chin Restraint is apparent, the Chin Restraint should be discarded and replaced.

CLEANING

When dirty, handwash in warm water using a mild detergent. Rinse thoroughly and drip dry.

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© 2018 ResMed Ltd. 168154/2 2018 02
Velcro is a registered trademark of Velcro Industries B.V.
Manufacturer: ResMed Ltd 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia