## ENGLISH

The information provided in this addendum replaces or supplements content found in the Stellar 100 and Stellar 150 User and Clinical Guides.

## Supplementary information on FiO<sub>2</sub> monitoring

**Note:** It may take up to 15 minutes for the FiO<sub>2</sub> sensor readings to meet the specified accuracy after powering on the device.

The FiO<sub>2</sub> sensor performance can be adversely affected by relative humidity, condensation on the sensor or unknown gas mixtures.

FiO<sub>2</sub> monitoring is affected by changes in ambient and cyclical pressure. This change in monitoring is proportional ie, at 3 000m (10 000 ft) altitude the ambient pressure decreases by 30% relative to sea level and will result in a 30% decrease in monitored FiO<sub>2</sub>.

For information on maintenance and disposal of the oxygen sensor, refer to the FiO<sub>2</sub> Monitoring Kit Guide for Stellar.

### FiO<sub>2</sub> sensor calibration

**Note:** Calibration must be done in room air without supplemental oxygen connection.

## Optional accessories

### Using water traps

When using a water trap, check the air circuit regularly for accumulated water. Follow the manufacturer's recommendations for frequency of inspection of water traps.

## Supplementary information on SpO<sub>2</sub>

Both the SpO<sub>2</sub> and pulse rate are logged at a sample rate of 1 Hz.

The optional SpO<sub>2</sub> sensor is calibrated for the display of functional oxygen saturation.

The Disconnected XPOD alarm will activate if the pulse oximeter is connected during therapy and an interruption to its data transfer occurs for more than 10 seconds.

For part numbers of compatible oximeter accessories for this device, see the ventilation accessories on [www.resmed.com](http://www.resmed.com) under the Products page. For information on how to use these accessories, refer to the user guide that comes with these accessories.

### WARNING

- Pulse oximeter sensors must not be used with excessive pressure for prolonged periods as this can cause patient pressure injury.
- The pulse oximeter sensor and cable needs to be verified for compatibility with Stellar, otherwise, patient injury can result.

## Internal battery

The internal battery will operate for approximately three hours under normal conditions.

It can take up to four hours to fully recharge the internal battery, however, this can vary depending on environmental conditions and if the device is in use.

## Working with alarms

When an alarm is muted, the Alarm mute key LED will light constantly. If after two minutes the problem is still present, the alarm will sound again for high or medium priority alarms except for the **Internal Battery Low (30%)** alarm. This alarm will be permanently muted with both the audible signal and alarm message cleared. However, these indicators will initiate again when the **Internal Battery Empty (2 minutes remaining)** alarm is activated.

## Technical specifications

Alarm Volume Range	High, medium and low priority alarm (measured according to IEC 60601-1-8: 2012): minimum 42 dB maximum 79 dB
Applied parts (as per IEC 60601-1)	Patient interface (Mask, endotracheal tube or tracheostomy tube) Oximeter (optional accessory)
FiO <sub>2</sub> measurement	Range: 18% to 100%, Resolution: 1%, Accuracy: $\pm(2.5 \text{ vol.}\% + 2.5\% \text{ of actual oxygen concentration})^*$ (as per ISO 80601-2-55)
SpO <sub>2</sub> measurement	Range: 70% to 100%, Resolution: 1%, Accuracy: $\pm 2\%$ in a range of 70% to 100%
Internal battery	Lithium-Ion battery, 14.4 V, 2.75 Ah, 40 Wh Operating hours: 3 h with a new battery under normal conditions

\* The measured FiO<sub>2</sub> sampled at a minimum of 1Hz has a response time of 14s to 90% of the final value. The long term output drift in 100% O<sub>2</sub> is typically <5% over 1 year.

## General warnings and cautions

### WARNING

Single-use products are not to be re-used due to degradation and contamination risks.

## Symbols

 Device weight;  Contains no China environmental hazardous substances.

## Guidance and manufacturer's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment –guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles. Single phase: at 0° 0 % UT; 250/300 cycle	0% UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles. Single phase: at 0° 0 % UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power source.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment –guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz  6 Vrms in ISM bands between 150 kHz and 80 MHz	3 Vrms 150 kHz to 80 MHz  6 Vrms in ISM bands between 150 kHz and 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  Recommended separation distance:  $d = 0.35 \sqrt{P}$ outside ISM bands $d = 1.2 \sqrt{P}$ inside ISM bands $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz  where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup>  Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz (Professional Healthcare)  10 V/m 80 MHz to 2.7 GHz (Home Healthcare)	3 V/m 80 MHz to 2.5 GHz (Professional Healthcare)  10 V/m 80 MHz to 2.7 GHz (Home Healthcare)	

Notes:

- Ut is the AC mains voltage prior to application of the test level.
  - At 80 MHz and 800 MHz, the higher frequency range applies.
  - These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

## Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz to 80 MHz Outside ISM bands	150 kHz to 80 MHz Inside ISM bands	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = 0.35 \sqrt{P}$	$d = 1.2 \sqrt{P}$	$d = 1.2 \sqrt{P}$	$d = 2.3 \sqrt{P}$
0.01	0.1	0.1	0.1	0.2
0.1	0.1	0.4	0.4	0.7
1	0.4	1.2	1.2	2.3
10	1.1	3.8	3.8	7.3
100	3.5	12.0	12.0	23.0

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note:

- At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### Potential impact of electromagnetic disturbances

The loss or degradation of the following clinical functions due to electromagnetic disturbances could result in compromised patient safety:

- Accuracy of ventilation control
- Accuracy of monitoring of airway pressure and volume
- Therapy alarms

Detection of this degradation could be observed by the following device behaviour:

- Erratic ventilation delivery
- Rapid fluctuations in monitored parameters
- False activation of therapy or technical alarms

### WARNING

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

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For patent and other intellectual property information, see ResMed.com/ip. © 2019 ResMed Ltd. 248950/1 2019-06
See ResMed.com for other ResMed locations worldwide. Stellar is a trademark and/or registered trademark of the ResMed family of companies.
**MANUFACTURER** **ResMed Ltd**