

Humidifier/device compatibility list

This supplement specifies the devices compatible with ResMed's humidifiers. For full details regarding the correct use of these devices, please refer to the relevant User Guides.

Humidifier Device	Passive Humidifiers		Heated Humidifiers				
	Passover	HumidAire 2iC™	HumidAire 2i™	HumidAire 3i™ H3i Plus™	H4i™	H4i™ for Stellar	H5i™
AutoSet CS™ 2		✓	✓				
AutoSet™ CS (S9)							✓
AutoSet Spirit™	✓	✓	✓				
Auto 25 (S9)							✓
S7™ Series	✓	✓	✓				
S8 AutoScore™							
S8 AutoSet Vantage™				✓	✓		
S8 AutoSet Spirit™							
S8 Escape™	✓			✓	✓		
S8 Lightweight™							
S8 Escape™ II				✓	✓		
S8 Elite™ II							
S8 AutoSet™ II				✓	✓		
S8 AutoSet Spirit™ II							
S9 AutoSet™ / S9 Elite™							✓
S9 Escape™ / S9 Escape™ Auto							
VPAP™ III Series	✓	✓	✓				
VPAP™ III ST-A with QuickNav		✓	✓				
VPAP™ Adapt SV		✓	✓				
VPAP™ Adapt (S9)							✓
VPAP™ Auto				✓	✓		
VPAP™ Auto (S9)							✓
VPAP™ S (S9)							
VPAP™ ST (S9)							
VPAP™ Auto 25							
VPAP™ ST					✓		
VPAP™ IV							
VPAP™ Tx			✓				
VPAP Malibu™			✓				
Stellar™						✓	

Note: Not all devices are available in all regions.