



Mirage™ SoftGel

NASAL MASK

Instant Softness, Lasting Comfort

The Mirage SoftGel Nasal Mask offers unprecedented softness and comfort. The DoubleGel cushion, with its clear gel layer provides comfort, while the blue gel layer maintains support and stability.

Backed by ResMed's renowned performance and quality, the frame accommodates two interchangeable cushions—Mirage SoftGel and Mirage Activa™ LT.

- **Soft, inviting, makes a great 'first impression'**
- **Excellent for new users**
- **Comfort throughout therapy**
- **Renowned quality and reliability**

FOR THE PATIENT

Soft comfort meets stability

- DoubleGel cushion combines softness with stability

A good night's sleep

- An array of small, built-in vents diffuses air quietly providing a peaceful environment for patients and their bed partners.

Fitting made easier

- The Mirage dual-wall with the DoubleGel cushion makes fitting simple.
- Cushion conforms to the unique facial contours of each individual, creating a comfortable fit.

Interchangeable cushions offer choice

- Choice of cushions makes it easier to remain compliant without completely replacing the mask.

FOR THE PROFESSIONAL

Customer satisfaction, added compliance

- Increased therapy acceptance through first-time fit—lasting comfort.
- Comprehensive user support package.

Improved efficiencies

- Interchangeable cushions and shared components improve inventory efficiencies, reducing costs, fitting and training times.
- A multi-patient, multi-use mask, Mirage SoftGel is suitable for repeated use in the lab.
- Extremely durable, Mirage SoftGel ensures consistent performance with fewer callbacks.

PATIENT-PREFERRED

- 8 out of 10 patients rate Mirage SoftGel better in comfort and durability than other gel masks*.
- 3 out of 4 patients rate their sleep better with Mirage SoftGel*.

* ResMed User Evaluation Trials—Mirage SoftGel vs Competitor Gel Masks



RESMED

MicroFit dial

offers 24 positions for better fit, seal and comfort

ResMed's unique frame

accommodates both Mirage SoftGel and Mirage Activa LT cushions for business efficiency, cost savings and flexibility

360° elbow rotation

positions the tubing for ease of use


Quiet, diffused vent design

disperses air softly and quietly, away from bed partner



Elbow partition

aids CO₂ washout and reduces noise output



Mirage SoftGel's unique DoubleGel cushion

Clear gel layer provides softness and comfort; blue gel layer provides support and stability

Set-and-forget headgear clips

retain optimal headgear settings; easy to attach and detach

Swivel

offers a convenient way to connect and disconnect air tubing

PRODUCT CODES

Mirage SoftGel Nasal Mask

Europe ¹	Europe ²	Europe ³
61605 (S)	61636 (S)	61610 (S)
61606 (M)	61637 (M)	61611 (M)
61607 (L)	61638 (L)	61613 (L)
61608 (LW)	61639 (LW)	61614 (LW)

¹ English (UK), French, German and Italian (Switzerland)

² English, Danish, Finnish, Norwegian and Swedish

³ English, Czech, Dutch (Belgium), French (Belgium), German (Belgium), Greek, Hungarian, Italian, Polish, Portuguese, Russian, Slovak, Spanish and Turkish

M Medium

L Large

St Standard

LW Large Wide

Manufacturer: ResMed Ltd 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia. Distributed by: ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA; ResMed (UK) Ltd 96 Milton Park Abingdon Oxfordshire OX14 4RY UK. See www.resmed.com for other ResMed locations world-wide. Activa and Mirage are trademarks of ResMed Ltd and Mirage and Mirage Activa are registered in U.S. Patent and Trademark Office. © 2009 ResMed Ltd. Specifications may change without notice. 1013843/1 09 11