

# RESMED

For the ultimate in stylish contemporary design, the S8 system immediately reflects its heritage and quality. When integrated with a ResMed mask and HumidAire 3i™, the S8 is the smallest truly integrated system available—so it sits attractively (and unobtrusively) by your bed as you sleep peacefully. In fact, the S8 is so compact, quiet and reliable, you will hardly know it's there.

## S8™ System



## S8™ System

Patient Information

**Unobtrusive in your  
bedroom—unobtrusive  
in your life**



ResMed Corp Poway, CA, USA +1 858 746 2400 OR 1 800 424 0737 (toll free), ResMed Ltd Bella Vista, NSW, Australia, +61 (2) 8884 1000 or 1 800 658 189 (toll free). Offices in Austria, Finland, France, Germany, Hong Kong, Japan, Malaysia, Netherlands, New Zealand, Singapore, Spain, Sweden, Switzerland, United Kingdom (see website for details). US Patents: 4944310 5199424 5245995 5522382 5704345 6029665 6138675 6240921 6363933 6367474 6398739 6425395 6502572 6591834 6635021 6705315 6745768. Others pending. Design registrations pending. S8, HumidAire 3i and Mirage Swift are trademarks of ResMed Ltd and HumidAire and Mirage are registered in the US Patent and Trademark Office. Specifications may change without notice. ©2005 ResMed Ltd. 1010342/1 05 06



## ResMed's S8™ System

### Compact without compromise

Small yet perfectly formed, ResMed's S8™ System redefines the shape of OSA therapy.

The system gives you uncompromising therapy in a compact, robust design—all in the palm of your hand.

### Humidification—it doesn't get easier than this

The HumidAire 3i offers unrivalled humidification in a compact device.

Its stylish design complements the contemporary look and feel of the S8, while its fliptop lid design and easy setup provides a therapy that is so easy to use.

Each model in the S8 system can be combined with the HumidAire 3i humidification module.

### Why humidify?

The cool, dry air from therapy may make your nose, mouth or throat dry, particularly in winter. It may also cause a runny nose, sneezing or nasal congestion. Heated humidification can help relieve these symptoms by adding moisture to the air before you breathe it. As well as improving comfort, humidification can help to prevent mouth breathing. Mouth breathing may decrease the effectiveness of your therapy, so it is important to manage this effectively.



### All packed and ready to go ...

The S8 brings you the freedom and flexibility you've been dreaming of.

Designed to be the most compact yet robust system on the market, it is incredibly convenient for travel. No matter where you are—from the beach to the peaks—ResMed's S8 has everything you need in a fully integrated package:

- exclusive tailored carry bag for flexible travel options
- altitude compensation for more accurate pressure delivery
- universal mains power supply automatically switches to local power
- 12V or 24 V DC input via ResMed's DC-12 converter.



### The complete therapy solution from ResMed

Having a properly fitted mask that is comfortable and doesn't leak is vital for successful OSA therapy. The S8 is fully compatible with the complete range of ResMed masks (except non-vented). Choose ResMed's Mirage Swift Nasal Pillows™ for one of the smallest integrated systems available.

With a reputation for comfort and seal, ResMed masks are available in a range of models and sizes to suit all face shapes.

