



FAA Compliance for ResMed Devices

To whom it may concern,

The US Department of Transportation (DOT) Final Rule, “*Nondiscrimination on the Basis of Disability in Air Travel*” (73 FR 27614 which updates Title 14 CFR Part 382), effective May 13, 2009 provides important new requirements for the accommodation of passengers with respiratory assistive devices (Ventilators, Respirators and CPAP machines).

In line with these requirements, respiratory assistive devices may be used onboard an aircraft, without further testing by the carrier, provided they have been tested for Electromagnetic Compatibility (EMC) in accordance with the current version of RTCA/DO-160, Section 21, Category M.

Some airlines may require notification greater than 48 hours before travel, and devices may need to be operated by battery. ResMed recommends that customers check with their airline.

ResMed has now successfully completed testing for the following respiratory assistive devices. **The devices listed below comply with RTCA/DO-160G, Section 21, Category M, and are considered FAA compliant.**

Product Code	Product Description
33030	S8 Compact
33001, 33007, 33013, 33036, 33060	S8 Escape
33032	S8 Lightweight
33033	S8 AutoScore
33048, 33051, 33052, 33053, 33054, 33061	S8 Escape II
33055	S8 Lightweight II
33021	S8 Elite
33039, 33045, 33062, 33072, 33073, 33074, 33075	S8 Elite II
33044	S8 AutoScore II
33112	S8 AutoSet Vantage
33105	S8 AutoSet
33128, 33137	S8 Respond
33129, 33141, 33150	S8 AutoSet II
33136, 33143, 33144, 33145, 33146, 33147, 33148	S8 AutoSet Spirit II
33138	S8 AutoSet C
26101, 26121	VPAP Auto 25
26119, 26120	VPAP S
26110, 26122	VPAP ST
26104, 26105, 26125, 26126	S8 Auto 25
26102, 26103, 26106, 26107, 26108, 26109, 26123, 26127	VPAP IV
26112, 26113, 26114, 26115, 26116, 26117, 26118, 26124	VPAP IV ST
36003, 36013, 36023, 36103, 36113, 36123, 36143, 36203, 36223, 36243, 36263, 36303, 36343, 36363	S9 Elite Series
36005, 36015, 36025, 36105, 36115, 36125, 36145, 36205, 36225, 36245, 36265, 36305, 36325, 36345, 36365	S9 AutoSet Series
36001, 36011, 36021, 36141, 36201, 36221, 36261, 36301, 36361	S9 Escape
36002, 36012, 36022, 36302, 36362	S9 Escape Auto
36004	S9 VPAP S - AMER
36006	S9 VPAP AUTO - AMER
36007	S9 VPAP ADAPT - AMER

36008	S9 VPAP ST - AMER
36014	S9 VPAP S + H5i - AMER
36016	S9 VPAP AUTO + H5i - AMER
36017	S9 VPAP ADAPT + H5i - AMER
36018	S9 VPAP ST + H5i - AMER
36024	S9 VPAP S + CLIM CNTRL - AMER
36026	S9 VPAP AUTO + CLIM CNTRL - AMER
36027	S9 VPAP ADAPT + CLIM CNTRL - AMER
36028	S9 VPAP ST + CLIM CNTRL - AMER
36100	S9 AUTOSET CS - FRA
36284	S9 VPAP S - FRA
36106	S9 AUTOSET 25 - FRA
36108	S9 VPAP ST - FRA
36110	S9 AUTOSET CS - TAKE AIR - DEU
36114	S9 VPAP S - TAKE AIR - DEU
36116	S9 AUTOSET 25 - TAKE AIR - DEU
36118	S9 VPAP ST 22 - TAKE AIR - DEU
36120	S9 AUTOSET CS - DEU
36124	S9 VPAP S - DEU
36126	S9 AUTOSET 25 - DEU
36128	S9 VPAP ST 22 - DEU
36140	S9 AUTOSET CS - GBR
36144	S9 VPAP S - GBR
36146	S9 AUTOSET 25 - GBR
36148	S9 VPAP ST - GBR
36200	S9 AUTOSET CS - EUR1
36204	S9 VPAP S - EUR1
36206	S9 AUTOSET 25 - EUR1
36208	S9 VPAP ST - EUR1
36220	S9 AUTOSET CS - EUR2
36224	S9 VPAP S - EUR2
36226	S9 AUTOSET 25 - EUR2
36228	S9 VPAP ST - EUR2
36304	S9 VPAP S - APAC
36360	S9 AUTOSET CS - APAC
36366	S9 AUTOSET 25 - APAC
36367	S9 VPAP ADAPT - APAC
36368	S9 VPAP ST - APAC
24921, 24923, 24925	ResMed Power Station II (RPSII)



Simon Lewi

Director – Regulatory Affairs, ResMed Ltd

22 November, 2011