



Mirage Kidsta™ Nasal Mask



The reassuring mask of choice for sleep-disordered breathing (SDB) in childhood

ResMed is the world leader in mask technology for SDB and we are bringing that expertise to a mask designed for children seven years and older.

The Mirage Kidsta™ Nasal Mask is based on the revolutionary range of adult Mirage™ designs. It combines a number of proven, practical and popular benefits with some child-friendly modifications. This ResMed mask is designed specifically for children who need continuous positive airway pressure (CPAP) or bilevel therapy.

ResMed's Kidsta Nasal Mask—with the proven and trusted Mirage dual-wall cushion technology—is a premium mask that combines comfort and performance in an unobtrusive design.



Stabilizing panel

holds headgear in correct position for optimal seal

Headgear

provides set and forget convenience

Mask elbow

connects tubing to mask frame. Provides greater ease for most children to click-on and click-off

Angle selector

enables customized fit for optimal comfort

Easy locking mechanism

provides easier handling and enables headgear clip to easily lock into mask frame

Technology highlight—leaders in mask technology

ResMed is renowned for its pioneering technologies and commitment to producing world-class products.

The Mirage Kidsta Nasal Mask was designed and manufactured based on anthropometric (facial and head) measurements on a range of children.

Data was assimilated and combined with ResMed's unique Mirage dual-wall cushion technology to help develop a mask cushion and headgear specifically designed for a child.



Benefits

- Dual-wall cushion reduces pressure on small noses and optimizes seal and comfort.
- Easy click-on, click-off parts are ideal for children's hands to grasp.
- Unique instruction booklet, specifically designed to support parent and child.
- Angle selector and stabilizing panels provide customized fit, enhanced stability, seal and comfort.
- Unobtrusive, streamlined design allows eye contact and clear vision. The open design also reduces feelings of anxiety and claustrophobia.
- Convenient set and forget headgear minimizes the need for continuous adjustment and allows greater control for smaller hands.
- Child can wear mask without elbow attached to tubing to facilitate acclimatization to therapy.



Open design provides child with a clear field of vision and encourages acclimatization.



The Mirage Kidsta Nasal Mask comes with a child's instruction booklet introduced by the mask mascot Oke. The guide includes a diary to help chart progress and is valuable when it comes to acclimatization and compliance.



Child-friendly design allows smaller hands to grasp parts easily.

Mirage Kidsta Nasal Mask System

Product Code 61011