

RESMED

Swift LT

NASAL PILLOWS SYSTEM

LESS IS MORE

- Minimal contact and optimal stability
- Freedom to move and sleep on your side
- Super-quiet and personalised comfort

Mask performance you can trust

ResMed's Swift masks have helped millions of users around the world discover the benefits of a nasal pillows system.

With Swift LT, less is certainly more. Its minimalist design masks many innovative features to enhance your experience of lightness and freedom—with stability.

Your mask comes with a comprehensive user support package, which includes a handy bag for storage and travel and a user guide with fitting and care instructions.



Swift LT—the latest generation nasal pillows from ResMed, leaders in mask technology

RESMED

Swift™ LT

NASAL PILLOWS SYSTEM

User Information



LESS IS MORE

ResMed Corp Poway, CA, USA +1 858 746 2400 or 1 800 424 0737 (toll free), ResMed Ltd Bella Vista, NSW, Australia +61 (2) 8884 1000 or 1 800 658 189 (toll free). (For the complete list of global offices, please see the ResMed website). Protected by patents: AU 785376, DE 20321468.4, DE 20321469.2, DE 20321470.6, DE 20321471.4, DE 20321472.2, EP 1314445, NZ 553756, NZ 553757, NZ 553822, NZ 553824, NZ 553825, US 6581594, US 6823865, US 7159587, US 7318437. Other patents pending. Protected by design registrations: AU 301943, CN 200430102469.3, JP 1243026, NZ 407161, US D552731. Others designs pending. Swift is a trademark of ResMed Ltd and Swift is registered in U.S. Patent and Trademark Office. © 2009 ResMed Ltd. 1012554/2 09 02

Minimal contact and optimal stability

- **Sheer minimalism**—Swift LT is unique among nasal pillows systems. Inspired by user feedback, 'less is more' goes beyond design—it's a total experience.
- **Minimal contact, optimal seal**—the first nasal pillows with ResMed's proven dual-wall technology. The pillows rest on a flexible cushion base maintaining seal, while you barely feel them on your face.
- **Light all the way through**—the entire mask provides a sense of 'less is more'. The headgear is shaped to feel secure, not tight. Lighter than standard tubes, the Swift LT tube minimises 'pull' on the mask.

Freedom to move and sleep on your side

- **Perfect for side-sleeping**—the streamlined narrow mask, made of ultra-soft material is unobtrusive when you sleep on your side.
- **Perfect for turning over**—headgear with flexible stabilising elements gently holds the system in place when you turn over on your side or on your back, and move during sleep. The new slim buckle at the back of the head is soft and designed for your comfort.
- **Keeps pace with your movement**—the cushion base enables the dual-wall pillows to move in different directions without compromising seal or creating leaks. This gives you more freedom to move.

Super-quiet and personalised comfort

- **Incredibly quiet**—this mask could well be the quietest* among nasal pillows, dispersing air gently through the fine vents. The innovative vent design directs air away from you and your partner to give you both a restful sleep.
- **Choice of tube direction**—the advanced design of the optional tube retainer enables you to choose the most convenient tube direction and keeps the tube in place there.
- **Click-fit comfort**—the click-fit adjuster selects the most comfortable angle for the nasal pillows to suit your profile. It adjusts the pillows in small increments until the mask feels just right for you.

* Only 25 decibels when tested as per ISO 3744:1994 at 10 cm H₂O

Swift™ LT

NASAL PILLOWS SYSTEM

LESS IS MORE

