


## Standard Trolley

Suitable for home and hospital use

### Purpose-built for ResMed's world-class range of ventilators

- Designed to carry ResMed's comprehensive range of ventilation devices—an interchangeable top plate to suit each model
- Perfect for the bedside—a basket neatly holds everything you need and an optional support arm keeps the air tubing away from the floor
- Stable and portable—with a stable five-wheel base (three are lockable) and an inbuilt shelf for the ResMed Power Station, the entire system can be transported on a single trolley
- Convenient humidifier bracket—makes it easy to hook on an external humidifier in one simple step



Some accessories are so well-designed, they become essential.



## USEFUL ADD-ONS

- Circuit support arm keeps the air tubing away from the floor, reducing any 'pull' on the patient's mask or tracheostomy
- Tube clips securely hold the air tubing in place
- Interchangeable top plate is available in a variety of sizes to suit various devices

**Premium circuit support arm**  
24917

**Top plate**  
For VS and Elisée™ Series  
24914  
For VPAP™ and Stellar™ Series  
24913

## PRODUCT CODES

|  |              |
|--|--------------|
| Standard trolley (without top plate)     | <b>24905</b> |
| VS and Elisée Series top plate           | <b>24914</b> |
| VPAP Series and Stellar Series top plate | <b>24913</b> |
| Premium circuit support arm              | <b>24917</b> |
| Tube clip—double                         | <b>24919</b> |
| Tube clip—single                         | <b>24920</b> |

**The new ventilator trolley from ResMed—it's more than an accessory; it's a necessary convenience.**

Note: Actual product available may vary from that shown in the picture.

**For other ResMed accessories and spare parts please visit [www.resmed.com/ventilationaccessories](http://www.resmed.com/ventilationaccessories).**