



# Data management/Device compatibility list

This supplement specifies the devices compatible with ResMed data management products. For full details regarding the correct use of these devices, see the relevant User Guides.

Note: Not all products are available in all regions.

	ResScan™ 6.0	S9™ Wireless†	EasyCare Tx™ / Tx Link™ 9.0	AirView™	RCM	RCMH**
<b>AirMini™</b>				✓		
<b>AirSense 10 / AirCurve 10 Series</b>						
AirSense™ 10 AutoSet	✓		✓	✓		
AirSense™ 10 AutoSet for Her	✓		✓	✓		
AirSense™ 10 Elite	✓			✓		
AirSense™ 10 CPAP	✓			✓		
AirCurve™ 10 S	✓			✓		
AirCurve™ 10 ST	✓			✓		
AirCurve™ 10 VAuto	✓			✓		
AirCurve™ 10 ASV	✓		✓	✓		
AirCurve™ 10 ST-A	✓			✓		
<b>S9 Series</b>						
S9 AutoSet™	✓	✓	✓	✓		
S9 AutoSet™ for Her	✓	✓		✓		
S9 Elite™	✓	✓	✓	✓		
S9 Escape™	✓	✓	✓	✓		
S9 Escape Auto™	✓	✓		✓		
S9 VPAP™ Tx			✓			
VPAP™ S (S9)	✓	✓		✓		
VPAP™ ST (S9)	✓	✓	✓	✓		
VPAP™ Auto (S9)	✓	✓	✓	✓		
VPAP™ Adapt (S9)	✓	✓	✓	✓		
VPAP™ Adapt (S9 ASV)	✓	✓	✓	✓		
VPAP™ ST-A (S9 iVAPS)	✓	✓	✓	✓		
VPAP™ COPD (S9)*	✓	✓		✓		
Lumis™ Tx	✓					
<b>S8 Series</b>						
S8 AutoSet™ II	✓		✓			
S8 AutoSet Vantage™	✓		✓			
S8 Elite™	✓		✓			
S8 Elite™ II	✓		✓			
<b>S7 Series</b>						
S7 AutoSet Spirit™	✓					
S7 Elite™	✓					
VPAP™ III Series	✓		✓			
VPAP™ Auto	✓					
VPAP™ Auto 25†	✓					
VPAP™ III ST-A with QuickNav	✓		✓			
VPAP™ Adapt SV	✓		✓			



	ResScan™ 6.0	S9™ Wireless†	EasyCare Tx™ / Tx Link™ 9.0	AirView™	RCM	RCMH**
VPAP Malibu™	✓					
VPAP™ S†	✓					
VPAP™ ST†	✓		✓			
<b>Ventilator</b>						
Astral™ 100	✓			✓	✓	✓
Astral™ 150	✓			✓	✓	✓
Stellar™ 100	✓		✓	✓	✓	
Stellar™ 150	✓		✓	✓	✓	
VS III™	✓					
<b>Home sleep testing device</b>						
ApneaLink™ Air				✓		

RCM – ResMed Connectivity Module

RCMH – ResMed Connectivity Module Hospital

† U.S. and Canada only.

\* U.S. only.

\*\* For a full list of available data for use by an EMR system, refer to the list on [www.resmed.com/downloads/devices/RCMH](http://www.resmed.com/downloads/devices/RCMH).

 **MANUFACTURER** ResMed Ltd 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia

See ResMed.com for other ResMed locations worldwide. AirMini, AirCurve, AirSense, AirView, Astral, ApneaLink, Lumis, S8, S9, Stellar, AutoSet, S9 Escape, S9 Elite, S9 Escape Auto, S8 Elite, S8 AutoSet Vantage, ResScan, VPAP Tx, S8 AutoSet, S9 AutoSet, VS III, VPAP Auto, VPAP Adapt SV and VPAP Malibu are trademarks and/or registered trademarks of the ResMed family of companies. For patent and other intellectual information, see ResMed.com/ip. © 2018 ResMed Ltd. 1018020/13 2018-05

ResMed.com