






Power supply unit/device compatibility list

English

This supplement specifies the power supply units (PSU) compatible with ResMed devices. For full details regarding the correct use of these devices, please see the relevant device User Guides.


Notes:

- Not all products are available in all regions.
- The PSU part number is located in the label at the back of the unit.

Device	Power supply unit				
	370001	370006	369102	R270-7198	370005
					
	Original PSU for Air10 and Lumis devices.	Compact PSU for Air10 devices.	Packed with the S9™ series.	Original PSU for Astral.	Compact PSU for Astral
AirSense™ 10 / AirStart™ 10	✓	✓	–	–	–
AirCurve™ 10 VAuto / AirCurve™ 10 S / AirCurve™ 10 ST	✓	✓	–	–	–
AirCurve™ 10 ASV / AirCurve™ 10 ASVAuto	✓	✓	–	–	–
AirCurve™ 10 CS Pacewave™	✓	✓	–	–	–
Lumis™ S / Lumis™ ST	✓	✓	–	–	–
AirCurve™ CS-A	✓	–*	–	–	–
AirCurve™ ST-A	✓	–*	–	–	–
Lumis™ ST-A	✓	–*	–	–	–
S9™	–	–	✓	–	–
ResMed Power Station II	✓	–*	✓	–	–
Stellar™	–	–	–	–	–
Astral™	–	–	–	✓	✓

*These device have not been validated for use with this PSU. The PSU will deliver the same power to the unit and meet the power requirements of the device, however the Electromagnetic Compatibility in the alarmed variant configuration and during charging of the RPSII has not been tested for.

Sample label:



ResMed Pty Ltd
 1 Elizabeth Macarthur Drive
 Bella Vista NSW 2153 Australia

See ResMed.com for other ResMed locations worldwide. AirSense, AirCurve, AirStart, Astral, Lumis, Pacewave, S9, and Stellar are trademarks and/or registered trademarks of the ResMed family of companies. For patent and other intellectual property information, see ResMed.com/ip. © 2020 ResMed. 10111949/1 2020-08