


AirMini™ Solutions



Rewrite your patients' bedtime stories

ResMed's AirMini solutions deliver therapy efficacy in a compact and portable package that fosters engagement and compliance.¹⁻² Small in size, big on features, AirMini solutions provide patients with the same therapy you've come to expect from ResMed in a simple, travel-friendly package.



[ResMed.com/AirMini](https://www.ResMed.com/AirMini)

AirMini: effective compact therapy

The result of our dedication to clinical research, innovation and development, AirMini solutions provide your patients with a complete but compact system. Leveraging the clinically proven³⁻⁶ algorithms used in our popular bedside devices, AirMini delivers a complete solution that covers setup to long-term follow-up.

With embedded connectivity, AirMini solutions keep you and your patients connected to therapy and offer the potential to generate a range of time-saving efficiencies.



ActiveAir™

This innovative technology is designed to regulate airflow to keep pressure delivery stable and comfortable.

HumidX™ waterless humidification*

Compact and discreet, HumidX, HumidX Plus and HumidX F20 are good options for patients looking to experience therapy comfort on the go.



* HumidX and HumidX Plus are only compatible with AirFit N20, AirTouch N20, AirFit N30 and P10 for AirMini masks. HumidX F20 is only compatible with AirFit F20 and AirTouch F20 masks.



AirFit™ and AirTouch™ series of masks

AirMini is compatible with the N20 and F20 series of masks, as well as AirFit F30, P10 and N30 for AirMini*.

Start & Stop

With only one button, it's so easy to use the device every night.

Bluetooth connection

Just press the Bluetooth button to synchronise AirMini with a smart device.

AirMini App™ & AirView™

AirView and AirMini App will seamlessly keep your patient connected to therapy, from initial set up through early coaching and troubleshooting to long-term usage.



A system to meet your patients' wants and needs

The AirMini therapy system is very appealing to patients: it offers the benefits of a standard bedside system in a compact package. You can provide them with more options and benefits when and where it suits them – at home and on the go.

These options include various algorithms, therapy mode and features, as well as a full range of mask options (nasal pillows, nasal and full-face).



Boosting patient engagement

The AirMini App allows your patients to tailor their comfort settings and view their therapy progress easily.

Your patients can also allow you to follow up on their progress by sending their detailed therapy data to AirView.

The AirMini App helps patients adhere to therapy earlier and stay engaged. This limits the patient drop-out rate⁷, allowing you to acquire and service more patients while increasing the lifelong value of each patient to your business.



Accessories for mobile therapy

AirMini solutions include a range of accessories that promote independence and enable patients to easily continue therapy even when they're on the move. These include a mobility bag, a DC/DC converter and a mount system.

The DC/DC Converter 65W allows your patients to operate their machine from a vehicle's cigarette lighter socket.


The AirMini travel case offers protection, while a travel bag makes it possible to transport the entire AirMini system (tubing, mask and device) with ease.

The AirMini mount system allows the AirMini device to be secured in a variety of bedside configurations for greater convenience and discretion.

Technology that serves patients and you

To provide performance in such a small package, ResMed designed the innovative ActiveAir technology, the crucial pivot that makes AirMini solutions possible.

Because the venting technology is built into the tube, the AirMini device is only compatible with masks that are designed to work with this venting system: the AirFit F30, N20 and F20 series of masks, as well as the AirFit P10 and AirFit N30 for AirMini*.




An innovative way to deliver humidification

ResMed has developed a venting technology to deliver portable waterless humidification* and performance to your patients, all in a small, light and lifestyle-friendly package.

HumidX and HumidX Plus are only compatible with AirFit N30, AirFit N20, AirTouch N20 and AirFit P10 for AirMini masks. HumidX F20 is only compatible with AirFit F20 and AirTouch F20 masks.

* HumidX and HumidX Plus are only compatible with AirFit N30, AirFit N20, AirTouch N20 and AirFit P10 for AirMini masks. HumidX F20 is only compatible with AirFit F20 and AirTouch F20 masks.



- 1 AirMini Patient Preference Study. Data on file: ID B98963.
- 2 AirMini CPL report. Data on file: ID A4056017.
- 3 Netzel et al. APAP device technology and correlation with patient compliance. *Somnologie - Schlafforschung und Schlafmedizin* 2014; DOI: 10.1007/s11818-014-0662-0.
- 4 Isetta et al. Comparative assessment of several automatic CPAP devices responses: a bench test study. *ERJ Open Res* 2015; 1:00031-2015.
- 5 Zhu et al. All APAPs are not equivalent for the treatment of sleep disordered breathing: a bench evaluation of eleven commercially available devices. *J Clin Sleep Med* 2015 11(7):725-34.
- 6 Isetta V et al. Novel Approach to Simulate Sleep Apnea Patients for Evaluating Positive Pressure Therapy Devices. *PLoS One*. 2016;11(3):e0151530. Published 2016 Mar 15.
- 7 Price Waterhouse Coopers. Empowering the sleep apnoea patient. A study of how myAir users behave better. 2016.

See ResMed.com for other ResMed locations worldwide. AirFit, AirTouch, ActiveAir, AirMini and HumidX are trademarks and/or registered trademarks of the ResMed family of companies. For patent information, see ResMed.com/ip. © 2020 ResMed. 1019978/4 2020-10